

BETA-GALACTOSIDASE ENZYME ASSAY SYSTEM

This protocol is adapted from Promega's Beta-Galactosidase Enzyme Assay System Protocol by the Gene Expression Lab.

This protocol is for use with Promega's Beta-Galactosidase Enzyme Assay Systems. For additional technical inquiries, contact Technical Service at 800-356-9526 or www.promega.com

BEFORE STARTING THE EXPERIMENT β -GALACTOSIDASE ASSAY STANDARD CURVE

BEFORE STARTING THE EXPERIMENT

- Cells are collected from transfection plates following protocol in Transient Co-Transfection Protocol
 - Thaw 2X Buffer
 - Prepare 1X RLB (Reporter Lysis Buffer)
 - Add 4 volumes of water to 1 volume of 5X RLB to produce a 1X stock solution.
-

β -Galactosidase Assay

1. Place the Assay 2X Buffer on ice.
2. It will be necessary to dilute the cell lysates in 1X Reporter Lysis Buffer. A 2:1 dilution of lysate to 1X Reporter Lysis Buffer (100 μ l of lysate plus 50 μ l of 1X Reporter Lysis Buffer) is a good starting dilution, but up to 150 μ l of cell lysate can be used per reaction. As a negative control, prepare the same dilution of a cell lysate made from cells that have not been transfected with the β -galactosidase gene.
3. Pipet 150 μ l of the appropriately diluted (or undiluted) cell lysates into labeled tubes.
4. Add 300 μ L of DEPC H₂O into each tube (to bring the total volume for each tube of the assay to 1.1 mL for spectrophotometer requirements)
5. Add 150 μ L of Assay 2X Buffer to each of the tubes & Mix all samples by vortexing briefly.
6. Incubate the reactions at 37°C for 60 minutes or until a faint yellow color develops. Color development continues for approximately 3 hours. If enzyme activity is low, samples may be incubated overnight if the reaction tubes are tightly capped.

7. Stop the reactions by adding 500 μL of 1M Sodium Carbonate (1x RLB used to lyse cells) or 500 μL of 1m Tris Base (1x RLA used to lyse cells). Mix by vortexing briefly.
8. Read the absorbance at 420nm.

Standard Curves

1. If a standard curve is desired, use standards between 0 and 6.0×10^{-3} units of β -Galactosidase. Prepare the following dilution series in 1X Reporter Lysis Buffer immediately before use.
2. Add 10 μL of 1u/ μL β -Galactosidase to 990 μL of 1X Reporter Lysis Buffer and mix.
3. Then add 10 μL of this 1:100 dilution to 990 μL of 1X Reporter Lysis Buffer and mix it to make a 1:10,000 stock solution.
4. Using this stock, prepare 150 μL of each β -Galactosidase standard per tube as described below.

β -Galactosidase Standard (milliunits)	Volume of 1:10,000 Stock	Volume of 1X Reporter Lysis Buffer
0	0 μl	50 μl
1.0	10 μl	40 μl
2.0	20 μl	30 μl
3.0	30 μl	20 μl
4.0	40 μl	10 μl
5.0	50 μl	0 μl

5. Follow the protocol described in Steps 4–8 in the β -Galactosidase Assay section above.
6. Plot the absorbance at 420nm versus concentration of β -Galactosidase standards.
7. If a standard curve is used, prepare fresh enzyme dilutions.